



P.1.866.609.6827

F.1.866.614.5619

orders@lighthousemedicaltd.com

www.lighthousemedicaltd.com

## EASyS Advantage CAD Order Form 2.06.18

Date:	Contact:	PO#	Marked For:
Customer Name:			
Address:			
Confirmation Email To (Buyer):		Phone:	Fax:


### Step 1 - Seat Sizes (required options indicated by \* symbol)



	EASyS <input type="checkbox"/> Size 1 \$2095	EASyS <input type="checkbox"/> Size 2 \$2230	Description
Seat*	660961	670961	Solid contoured, integrated seat/back, 180 degree recline, growth adjustable; seat depth/ width adjustable from underneath, articulated foot plate height and angle adjustable; front/rear facing seat, easy to change; ventilated/washable upholstery; knee angle adjustment with release lever; cozy comfortable seat cushion, BSS Balance Safety System; calf strap; 5 different color options.
Back*	660962	670962	
Recline*	660963	670963	
Hip Sup.	660964	670964	
Calf Strap*	6635	6735	
Foot Plate*	660965	670965	

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Seat to Head	Max Height	Hip Angle
Size 1	7 - 11.5"	7.5 - 11.75"	7.5 - 11.75"	16.5 - 26.75"	49"	90 - 180 Degrees
Size 2	8.5 - 14.5"	10.5 - 15"	10 - 14.5"	21.5 - 29"	60"	90 - 180 Degrees

### Upholstery Color Options\*



Size 1	6247 <input type="checkbox"/>	6248 <input type="checkbox"/>	6249 <input type="checkbox"/>	6246 <input type="checkbox"/>	6245 <input type="checkbox"/>
Size 2	6447 <input type="checkbox"/>	6448 <input type="checkbox"/>	6449 <input type="checkbox"/>	6446 <input type="checkbox"/>	6445 <input type="checkbox"/>

### Step 2 - Frame Options (Required options are indicated below by \* symbol)

**Advantage-chassis** - Folding aluminum frame, tilt-in-space, spring suspension cargo basket.

	Size 1		Size 2	
Advantage-chassis	62023 <input type="checkbox"/>	\$2650	64023 <input type="checkbox"/>	\$2720
Castor Pin Lock*	67830 <input type="checkbox"/>	Included	67830 <input type="checkbox"/>	Included
Bus Tie-Down Loop Kit*	67839 <input type="checkbox"/>	Included	67839 <input type="checkbox"/>	Included
Pneumatic Wheels, Full Set	6646 <input type="checkbox"/>	\$260	6646 <input type="checkbox"/>	\$260
Frame Padding, Reversible	6622 <input type="checkbox"/>	\$155	6622 <input type="checkbox"/>	\$155
Umbrella with Advantage-chas	6594 <input type="checkbox"/>	\$115	6594 <input type="checkbox"/>	\$115
Tip Assist	6750 <input type="checkbox"/>	\$95	6750 <input type="checkbox"/>	\$95
Cargo Basket	Included		Included	



### Specifications

	Size 1	Size 2
Seat Capacity	77 lbs	88 lbs
Frame Capacity	99 lbs	110 lbs
Unit Weight	41 lbs	47 lbs
Seat Suspension	35 Degrees	
Upholstery	Trevira-CS	
Foot Angle	78 - 100 Degrees	
Footrest	11. - 7.012	12.5" x 8.25"
Floor to Seat	19"	19.5"
Max Width	23.25"	27.25"
Fold Length	32"	34"
Fold Height	16.5"	21.5"
Brakes	Foot Operated	
Seat Direction	Front and Rear Facing	
Front Wheels	7.5" solid	
Rear Wheels	10" solid	

**Q-chassis** - Hi-low indoor base/feeder/therapy frame, height/tilt adjustable, pedal controlled gas spring, lockable rear wheels

	Size 1		Size 2	
Q-chassis	6666 <input type="checkbox"/>	\$2295	6666 <input type="checkbox"/>	\$2295
EASyS Adapter for Q-chassis	62951 <input type="checkbox"/>	Included	62951 <input type="checkbox"/>	Included
Push Bar for Q-chassis	6660 <input type="checkbox"/>	\$370	6660 <input type="checkbox"/>	\$370



### Q-chassis Technical Specs\*\*\*

Capacity	88 LBS
Seat Tilt	-7 to 31 Degrees
Floor to Seat	11 to 28"
Seat Direction	Front Facing



## EASyS Advantage CAD Order Form 2.06.18

### Step 3 - Positioning and General Options

Lower Body Positioning	Unit	Small	Medium	Large	Price	Order
Foot Straps with Cam Buckle	Pair		6648 <input type="checkbox"/>			\$45
Ankle Huggers, Bodypoint	Pair		213151 <input type="checkbox"/>			\$250
Footrest Pad, Height 1 3/4	Each	10622 <input type="checkbox"/>		10023 <input type="checkbox"/>		\$115
<b>Mid Body Positioning</b>						
Abduction Block, Padded	Each	6620 <input type="checkbox"/>		6720 <input type="checkbox"/>		\$185
Hip Supports 2 Point, Padded SR Buckle	Each	8420 <input type="checkbox"/>		8470 <input type="checkbox"/>		\$90
Hip Supports 4 point, Padded SR Buckle	Each	8421 <input type="checkbox"/>		8471 <input type="checkbox"/>		\$245
Soft Pelvic Harness	Each	6614 <input type="checkbox"/>		6714 <input type="checkbox"/>		\$195
Crotch Strap (fits 6629, 8420, 8470, 8421, 8471)	Each		6341 <input type="checkbox"/>			\$30
Tray, Height Adj, Clear with Picture Window	Each	66281 <input type="checkbox"/>		67281 <input type="checkbox"/>		\$545
Mounting Hardware Tray * Req'd For 6628/6728	Each	66282 <input type="checkbox"/>		67282 <input type="checkbox"/>	Included	
<b>Upper Body Positioning</b>						
Chest Butterfly Style Neoprene	Each	6385 <input type="checkbox"/>		6585 <input type="checkbox"/>		\$215
Chest Vest Style	Each	6631 <input type="checkbox"/>	6632 <input type="checkbox"/>	6731 <input type="checkbox"/>		\$185
Chest H-Style, Posture Support, Bodypoint®	Each		214806 <input type="checkbox"/>			\$225
PivotFit Shoulder Harness, Dynamic, Bodypoint®	Each		214270 <input type="checkbox"/>			\$225
PivotFit Shoulder Harness, Bodypoint®	Each		214250 <input type="checkbox"/>			\$200
Trunk Support, Rigid Lateral	Pair	6634 <input type="checkbox"/>		6734 <input type="checkbox"/>		\$300
Trunk Support, Flexible Lateral with Chest Strap	Pair		6633 <input type="checkbox"/>			\$420
Pillow contoured	Each		6824 <input type="checkbox"/>			\$200
Head Support Rigid Laterals	Pair	6624 <input type="checkbox"/>		6724 <input type="checkbox"/>		\$225
<b>General Accessories</b>						
Grab Bar, Padded Black	Each	6621 <input type="checkbox"/>		6721 <input type="checkbox"/>		\$140
Vent Tray	Each		110004 <input type="checkbox"/>			\$525
IV Hanger With Mounting Clamp	Each		665011 <input type="checkbox"/>			\$200
Oxygen Tank Holder, E-Tank	Each		665012 <input type="checkbox"/>			\$200
Straps, Secure folded seat/frame	Each		667011 <input type="checkbox"/>			\$35
Straps, Wrist "Runaway" Tether	Each		667111 <input type="checkbox"/>			\$35
<b>Weather Protection</b>						
Canopy, Sun/Rain with Window Black	Each		663761 <input type="checkbox"/>			\$375
Lower Body Wind Breaker, black for 663761	Each		661852 <input type="checkbox"/>			\$130
Canopy, Rain Only, Clear	Each		6638 <input type="checkbox"/>			\$250
Rain Shield, Clear for Canopy 663761/637811	Each		6739 <input type="checkbox"/>			\$150
Rain cape, Grey	Each	6340 <input type="checkbox"/>		6540 <input type="checkbox"/>		\$200
Sleeping bag, Summer Weight, Black	Each	6619 <input type="checkbox"/>		6719 <input type="checkbox"/>		\$385
Sleeping Bag, Winter Weight, Black	Each	6618 <input type="checkbox"/>		6718 <input type="checkbox"/>		\$400
Umbrella, Angle Adjustable	Each		6594 <input type="checkbox"/>			\$115