

anchor

bathing chair



 **irclespecialty**



+1.866.609.6827



www.circlespecialty.com



info@lighthousemedicaltd.com



The Anchor provides a carefree bathing experience for children and caregivers & make safety and comfort a priority so your child looks forward to bath time.

We designed the Anchor knowing your child's needs are unique and constantly evolving. That's why you can make angle and height adjustments with one hand on either side of the chair, all while your child is seated. It is also depth adjustable allowing for growth and long-term use. These agile and individualized functionalities ensure longevity, comfort, and hassle free bath time **so every bubble counts.**



Angle Adjustable

back, seat, & calf support



Folds Flat

easy storage & transport

FEATURES

- Depth adjustable for growth
- Angle adjustable back, seat and calf support, can be configured for individual requirements
- Height adjustable
- Make one handed angle & height adjustments on either side of the chair
- Positioning straps and head support can be placed anywhere along length of the chair
- Folds flat for easy storage and transport

ACCESSORIES

01 Rolling Shower Base

AN2100

Adjustable shower base with swivel wheels and locking casters. Accommodates both chairs at multiple height, depth, and angle positions.



Product Dimension	47 L x 23.5 W x 15.5 H
Product Weight	9 lbs

02 Lateral Support

AN2110

Wedges and adjustable strap help maintain optimal trunk alignment.




03 XL Positioning Straps

AN2111

Set of XL trunk & pelvic straps.



ADJUSTMENTS



SPECIFICATIONS

	AN1200 X-Small	AN1400 Small	AN1600 Large
Seat Width	12"	14.5"	16"
Seat Depth	10.5" - 13"	10.5" - 13"	13"-16"
Back Height	18.5"	24"	32"
Back Width	14.5"	17"	18.5"
Calf Support Length	10.2"	11"	14"
Seat Height	5.75"-10.5"-14"-15"	5.75"-11"-15.5"-17"	5.75"-13"-18.25"-20.5"
Back Angle	0°, 30°, 60°, 90°	0°, 30°, 60°, 90°	0°, 30°, 60°, 90°
Calf Support Angle	30°, 60°, 90°	0°, 30°, 60°	0°, 30°, 60°
Product Weight	12.75 lbs	14.5 lbs	18 lbs.
Weight Capacity	155 lbs	155 lbs	155 lbs

NOTE:

We regretfully cannot accept returns for bathing and toileting products due to public hygiene concerns.