

Gel cushion

Unique adhesive surface!

FEATURES AND BENEFITS

The new ConfortMedic gel cushion was designed to increase comfort and safety on most bathing lift. It offers the following benefits for patients and caregivers:

Increase comfort

Increase safety


Unique adhesive surface

Does not float

Hypoallergenic and non-toxic

Easy to disinfect

Made of silicon



DIMENSIONS

16 x 16 5/8 x 1 in

WARRANTY

1 year against any manufacturing defect.

MAINTENANCE

Clean and disinfect after each use with non-corrosive disinfectant. Rinse cushion with water and dry it with a clean towel. Put the cushion on a clean and dry towel, bonding surface up. Please read the instruction notice included with the cushion.

MANUFACTURED BY: CONFORTMEDIC

1140 Neveu, Saint-Césaire, Québec J0L 1T0

Toll free: 1 800 442-2246 Fax: 450 469-0368

WWW.CONFORTMEDIC.CA