



Thomashilfen

ThevoSleepingStar



Bed time

Many children with special needs have sleep disorders or disturbances, affecting parents and other family members. These disorders include: issues with falling asleep and staying asleep, restlessness during sleep, as well as disorientation upon awakening. Sleep disorders affect a child's natural growth and development.



“I was at the end of my rope and burned out emotionally. Nothing really helped. I am so glad that thanks to her ThevoSleepingStar my daughter is more refreshed and relaxed.”



Improved sleep with the Sensory Support Mattress

- ✓ MiS Micro-Stimulation® supports mobility
- ✓ Improves body perception
- ✓ Reduces spasms and pain
- ✓ Calms the child for restorative sleep

Cover

removable and easy to clean for optimal hygiene

Wing suspensions

have a positive influence on the children's restless sleep by providing tiny movement impulses

Mattress

is especially tailored for children with special needs to provide comfortable sensation of lying.

Positioning suspension

soft sensation of lying → extended sleep phases


is optimally tailored to the sleep and lying requirements of children with special needs

Children with severe motor limitations **need additional movement impulses**

Severely immobile children or children in coma vigil need additional movement impulses.

A balanced ratio between relaxation and movement provides restorative sleep.


In addition, risk of pressure ulcers is minimized.

A photograph of a woman with dark hair hugging a young girl from behind. The woman is smiling and looking towards the camera. The girl is also smiling and looking towards the camera. They are both wearing light-colored clothing. The background is a soft, out-of-focus indoor setting.

“We came to appreciate the Active SleepingStar in the Children’s Hospital first. Now, we have it at home, as well. All of us are doing much better now.”


Improved sleep with the Sensory Support Active Mattress

- ✓ Reduces pressure
- ✓ Reduces pain and spasms
- ✓ Provides extra movement impulses
- ✓ Improves quality of sleep



Active movements


a simple menu selection offers a variety of movement programs. This provides stimulation individually adapted to the child, without disturbing noises



Wing suspensions

take in micro motions and return them to the child

It is your choice!



Shared features:

- sewn zipper
- washable at 65°C / 149°F
- suitable for tumble dryers
- autoclavable
- suitable for the re-use

features: Incontinence cover

- breathable
- wipes clean with disinfectant



Sizes: Thevo Therapy mattress for children


- ThevoSleepingStar • ThevoSleepingStar A

Thevo Therapy mattress for children are available from a width of
Crib: 27"x55"

Hospital: 36"x79"

Twin: 39"x75"

Questions about sizes? Contact us, we are happy to help!



Visit us on www.thevo.info!

Thomashilfen

Experts in restorative sleep for patients in any setting - healthcare facilities or at home

Thevo – Mobility Beds

Survey on the Effectiveness of the MCS Micro-Stimulation®-System Thevo-Adapt® "Sleeping Star"

MIS Internationale Fördergemeinschaft "Micro-Stimulation" e.V.

With the collaboration of:

Neurologie Spezialambulanz
Klinik für Neurologie, Universität Tübingen
Martin Behrens
Lehrstuhl für Geriatrie, Universität des Saarlandes
Christine Mühlhahn
Neurologische Universitätsklinik

Disability

Ways	Medical (%)	Physical and mental (%)
None	11	11
Only	11	11
Light	11	11
Medium	11	11
Severe	11	11
Very severe	11	11
Total	100	100

Disability severely affected in the type of disability

Ways	Female	Male	Total
None	11	11	22
Only	11	11	22
Light	11	11	22
Medium	11	11	22
Severe	11	11	22
Very severe	11	11	22
Total	100	100	200

Significance: χ^2 -test: $p = 0.000$, $p < 0.001$

in addition referred from various studies:
reference: χ^2 -test: $p = 0.000$, $p < 0.001$

Ways	Female	Male	Total
None	11	11	22
Only	11	11	22
Light	11	11	22
Medium	11	11	22
Severe	11	11	22
Very severe	11	11	22
Total	100	100	200

Brochures and studies are available for download at www.thevo.info.



Thomashilfen North America

309 South Cloverdale Street

Unit B 12

Seattle WA 98108, USA

P 206 763 0754

F 866.870.0801



Walkmühlenstraße 1 · D - 27432 Bremervörde · Germany

Phone: +49 (0)47 61 / 886-63 or -68

Fax: +49 (0)47 61 / 886-19 · www.thevo.info