Sit-u-Up

The Sit-u-Up elevates a person smoothly from lying down in bed to sitting up.





Technical Specifications

| Max user weight | 150 kg / 24 stone |
|---------------------------|--------------------|
| Min height | Flat |
| Total length | 75 cm / 29 ½" |
| Width | 70 cm / 27 ½ " |
| Total Product weight | 4.8 kg / 10 ½ lbs |
| Back/headrest length | 73 cm / 28 ¾" |
| Backrest angle adjustment | Flat to 65 degrees |
| Airflo Compressor | Mk3 |

Product code: MPCA130200

Easy to fit

- Sits on top of the mattress
- With no fixtures or fittings, it can be used on any size bed
- One airline connects the Airflo to the Pillowlift. Press one button to inflate

Improves wellbeing

- Stable and comfortable, a person can choose to be lifted on their back or side
- The Sit-u-Up sits on top of the mattress enabling it to be fitted to any bed and causes no disturbance to a second person in double bed
- Minimises moving and handling risks so reducing the risk of musculoskeletal injury for the carer
- Its portable and lightweight which means it can be rolled up for travelling

Safe to use

- Includes a structured backrest ideal for those who need extra support
- Stable and compact
- Quilt provides comfort and support
- Inflation and deflation can be stopped at any time, allowing freedom to choose the most comfortable position
- Anti-slip pillow holder

Simple to maintain

- Lightweight and portable, the Sit-u-Up separates into two parts
- For a routine clean, sponge clean with a detergent solution and dry thoroughly
- The covers can be machined washed
- Spare covers can be purchased

