

EASYS ADVANTAGE Pediatric Stroller

Please send to : orders@lighthousemedicaltd.com

Phone 1.866.609.6827 or FAX to 1.866.614.5619



Date:	Contact:	PO#	Marked For:
Customer Name:			
Address:			
Confirmation Email To (Buyer):		Phone:	Fax:

EaSyS Advantage paediatric stroller has a solid contoured, integrated seat/back with 180 degrees recline and growth adjustable seat depth. ELR with calf pad is standard and articulating footplate is height and angle adjustable. Front/Rear facing seat is easy to change. Upholstery is ventilated and washable and each stroller has a comfortable seat cushion. Knee angle adjustment is made easy with the release lever and each unit has the BSS Balance Safety System. Available in two sizes and five different colour options.

Seat Size Chart		Lower Leg	Seat Depth	Seat Width	Back Height #	Max Height	Hip Angle
	Size 1	7 - 11.5"	7.5 - 11.75"	7.5 - 11"	16.5 - 26.75"	49"	90 - 180 Degrees
	Size 2	8.5 - 14.5"	10.5 - 15"	10.5 - 14.5"	21.5 - 29"	60"	90 - 180 Degrees

Model and Size:	Seating assembly includes the following:	MSRP	
EaSyS Advantage Size 1	660961 Seat 660962 Back 660963 Recline 660964 Hip 660965 F/plate 660966 ELR 6635 Calf panel	\$2800	
EaSyS Advantage Size 2	670961 Seat 670962 Back 670963 Recline 670964 Hip 670965 F/plate 670966 ELR 6735 Calf panel	\$2800	

UPHOLSTERY COLOUR -	Size 1/2	MUST Choose One
Red/Grey	6247/6447	
Blue/Grey	6248/6448	
Green/Grey	6249/6449	
Magenta/Grey	6246/6446	
Grey/Grey	6245/6445	

Choose Chassis Type - A or Q

A-Chassis - Folding aluminum frame, tilt-in-space, dual spring suspension basket, tilt-in-space lever at push bar for easy access

A- Chassis Size 1	#6202	\$3190	
A- Chassis Size 2	#6402	\$3190	
Castor Pin lock one size	#67830	included	
Bus Tie down Loop Kit	#67839	included	
Cargo Basket		included	
Pneumatic Wheels, full set one size	#6646	\$300	
Frame Padding	#6622	\$165	
Umbrella	#6594	\$140	
Tip Assist	#6750	\$130	

Q-Chassis - Hi-low Indoor base/feeder/therapy frame, height/tilt adjust w/pedal gas spring tilt, lockable rear wheels	MSRP	
Q Chassis One Size	#6666	\$2990
Adapter for QChassis	#62951	included
Push bar for Q-chassis One Size	#6660	\$440

Q-chassis Technical Specs***

Capacity	88 LBS
Seat Tilt	-7 to +31°
Floor to Seat	11 to 28"
Seat Direction	Front Facing



Thomashilfen
North America

Positioning and General Options							
	Lower Body Positioning	Unit	Small	Medium	Large	MSRP	Order
	Foot Straps with Cam Buckle	Pair	6648			\$55	
	Ankle Huggers, Bodypoint®	Pair	213151			\$300	
	Footrest Pad, Height 1 3/4	Each	10622		10023	\$145	
Mid Body Positioning							
	Abduction Block, Padded	Each	6620		6720	\$225	
	Hip Supports 2 Point, Padded SR Buckle	Each	8420		8470	\$140	
	Hip Supports 4 point, Padded SR Buckle	Each	8421		8471	\$190	
	Hip Supports Soft Pelvic Harness	Each	6614		6714	\$225	
	Crotch Strap (fits 8420, 8470, 8421, 8471)	Each	6341			\$45	
	Tray, Height Adj, Clear with Picture Window	Each	66281		67281	\$595	
	Mounting Hardware Tray * Req'd For 6628/6728	Each	66282		67282	*included with above	
	Arm Rest, Height, Depth, Angle Adjustable	Pair	8431		8478	\$525	
Upper Body Positioning							
	Chest Butterfly Style Neoprene	Each	6385		6585	\$270	
	Chest Vest Style	Each	6631	6632	6731	\$270	
	Chest H-Style, Posture Support, Bodypoint®	Each	14813			\$285	
	PivotFit Shoulder Harness, Dynamic, Bodypoint®	Each	214270			\$285	
	PivotFit Shoulder Harness, Bodypoint®	Each	214250			\$270	
	Rigid Lateral, left & right, height adjustable	Pair	6634		6734	\$355	
	Flexible Lateral, left & right, w/chest strap adj.	Pair	6633			\$420	
	Head Support, contoured Pillow	Each	6824			\$255	
	Head Support rigid lateral, left & right, adjustable	Pair	6624		6724	\$305	
General Accessories							
	Grab Bar, Padded Black	Each	6621		6721	\$180	
	Footrest tube upholstery, grey	Pair	6230			\$40	
	Vent Tray	Each	6294		6494	\$595	
	Oxygen Tank Holder, E-Tank	Each	665012			\$285	
	Straps, Secure folded seat/frame	Each	667011			\$45	
	Straps, Wrist "Runaway" Tether	Each	667111			\$45	
Weather Protection							
	Canopy, Sun/Rain with Window Black	Each	663761			\$465	
	Lower Body Wind Breaker, black for 663761	Each	661852			\$140	
	Canopy, Rain Only, Clear	Each	6638			\$355	
	Rain Shield, Clear for Canopy 663761/637811	Each	6739			\$185	
	Rain cape, Grey	Each	6340		6540	\$240	
	Sleeping bag, Summer Weight, Black	Each	6619		6719	\$430	
	Sleeping Bag, Winter Weight, Black	Each	6618		6718	\$440	
	Umbrella, Angle Adjustable	Each	6594			\$140	

Technical Specs	Size 1	Size 2
Seat Capacity	77 lbs	88 lbs
Frame Capacity	99 lbs	110 lbs
Unit Weight	36 lbs	42 lbs
Max Seat tilt	31"	
Upholstery	Trevira-CS	
Foot Angle	78" - 100"	
Size of Foot Rest	11.5" x 7"	13" x 8"
Floor to Seat	19"	19.5"
Max Width	23"	27"
Fold Length	32"	35.5"
Fold Height	18.5"	18.5"
Brakes	Foot Operated	
Seat Direction	Front and Rear Facing	
Front Wheels	7.5" solid	
Rear Wheels	10" solid	