

# PROTEKT® DASH

## 32500 | STANDING TRANSFER

The Protekt® Dash is your best choice for fall prevention, reducing caregiver lifting injuries and promoting patient mobility. The Protekt® Dash offers ample weight capacity at 500 lbs. and includes an adjustable base which provides clearance around common obstacles such as recliners allowing access to patients.

- **Dual shin pads for extra support.**
- **Molded seat can be rotated upwards to allow for patient loading.**
- **Cross bar for extra support.**
- **Non-Slip footplate for safety.**
- **Ergonomically designed foot pedal offers hands free opening of base.**
- **Dual rear locking casters.**
- **Simple solution to compliment your electric stand assist lifters.**
- **500 lb. weight capacity.**
- **First choice for any fall prevention program.**
- **Support strap available - attachment hooks for optional side support and back strap.**

### SPECIFICATION FOR ITEM # 32500

Weight Capacity	500 lbs.
Unit Weight	78 lbs.
Minimum Base Width	26.8"
Maximum Base Width	37.0"
Open Base Width	37.0"
Closed Base Width	26.7"
Base Opens	Yes
Overall Height	43.3"
Base Height	4.5"
Minimum Seat Height	26.7"
Seat Width	20.0"
Overall Length	35.4"
Knee Pads/Shin Guards	Yes
Optional Support Strap	Yes
Optional Seat Locks	Yes

**WARRANTY:**  
Frame 3 years



Distributed in Canada by  
LHM - Lighthouse Medical Ltd.  
1.866.609.6827  
[www.lighthousemedicaltd.com](http://www.lighthousemedicaltd.com)