



MANGAR ELK Lifting Cushion

The ELK is a lifting cushion designed to lift from the floor.

EASY TO FIT

- Lifting cushion is inflated using an Airflo Compressor with a simple hand control
- Simple to use, requiring very little training
- Suitable for use inside or out
- Can be used by one carer

SIMPLE TO MAINTAIN

- Clean using a proprietary non-abrasive mild liquid cleaner.
- Easy to deflate, fold away and store
- The material is designed to exceptionally high standards with extraordinary durability, strength & performance in the toughest of environments

SAFE TO USE

- Lifts the frailest person up to 980 lbs from the floor
- Exceptionally rigid and stable
- Ideal for carers of individuals with neurological disorders who fall regularly

IMPROVES WELL-BEING

- Provides a safe, dignified lift for a fallen person
- Reduces the risk of musculoskeletal injury for the carer
- Can be used in confined spaces where traditional lifting equipment may not fit

TECHNICAL SPECIFICATIONS

Max user weight	450 kg
-----------------	--------

Max height	56 cm
------------	-------

Min height	Flat
------------	------

Seat depth	50 cm
------------	-------

Width	57"
-------	-----

Total product weight	3.6 kg
----------------------	--------

Airflo Compressor	24
-------------------	----

Lifting performance per charge	Up to 9 lifts
--------------------------------	---------------



Features

1. **4-Way Hand Control.**
2. **Made from very strong polyurethane** coated nylon material.
3. **Remove or keep stability bar** to provide extra comfort for the user.
4. **Can be used to lift a stretcher.**
5. **Easy to setup** - connect the four hoses from the hand control to the corresponding numbers and coloured connectors on the ELK.
6. **Powered by an Airflo 24 Compressor.**

Accessories

- **ELK Stretcher Bars**
- **Slidesheet**
- **ELK End Cap Set**



Product code: MPCA070400



For more information visit

Lighthouse Medical Ltd.
Email: info@lighthousemedicaltd.com
Phone: 604.472.0119

WIN  CARE