



Knee Trekker Product Line



# Stay Active + Mobile using a KNEE TREKKER Knee Walker

Knee Walkers are an excellent short and long term alternative to crutches for those patients with lower leg injuries.



## EXPLORER

If you spend time outdoors, an all-terrain walker may be the right choice for you.



## EXCAPE

Excellent choice for both indoor and outdoor use.



## EXCURSION

Upgraded Econo version for the budget-minded person.

# KNEE EXPLORER TREKKER



## FEATURES



Rear Drum brakes increase the longevity of the tire and requires less maintenance than traditional systems



Automotive style tie rod steering system offers maximum stability and improved safety




Large 12" soft, foam filled tires provide a smooth ride on all terrains and are maintenance free



2" Premium Contoured Knee Pad



Folds for easy transport



Knee stem is height adjustable for users 4'6" - 6'6"



Handle bars are height adjustable 33"-42"



Basket, LED Light and bell included

# KNEE EXCAPE TREKKER



OPTIONAL WALKER BAG  
\*NOT INCLUDED

## FEATURES



Rear Drum brakes increase the longevity of the tire and requires less maintenance than traditional systems



Automotive style tie rod steering system offers maximum stability and improved safety



8" Solid Rubber tires are non-marking and maintenance free



2" Premium Contoured Knee Pad



Folds for easy transport



Knee stem is height adjustable for users 4'6" - 6'6"



Handle bars are height adjustable 33" - 42"



Basket, LED Light and bell included

# KNEE EXCURSION TREKKER



## FEATURES



Rear Drum brakes increase the longevity of the tire and requires less maintenance than tradition systems



8" Solid Rubber tires are non-marking and maintenance free



2" Premium Contoured Knee Pad



Folds for easy transport  
(Some disassembly required)



Knee stem is height adjustable for users 4'6" - 6'6"



Handle bars are height adjustable 33"-42"

# Compare Models

KNEE WALKERS	EXPLORER	EXCAPE	EXCURSION
GENERAL	9277	9276	9270
Weight (lbs)	23	22	21
Weight Capacity (lbs)	330	300	300
Recommended User Height	4'6" to 6'6"	4'6" to 6'6"	4'6" to 6'6"
Steerable	Yes	Yes	Yes
Folds For Easy Transport	Yes	Yes	Yes - Dissembly Required
Colour	British Racing Green	Golden Champagne	Glossy Black
ALL DIMENSIONS			
Overall Length	35"	30"	31"
Width	20"	19"	16"
Front Anxle Width	20"	19"	16"
KNEE PLATFORM			
Lowest Position	18"	17"	16"
Highest Position	24"	22"	22"
KNEE PAD			
Width	7"	7"	7"
Length	14"	14"	12.5"
Thickness	2.5"	2.5"	2"
HANDLE BARS			
Min Height	35"	33"	33"
Max Height	45"	42"	42"
WHEELS			
Diameter	12"	8"	8"
Number	4	4	4
Type	Solid Foam	Solid Non-Marking Rubber	Solid Non-Marking Rubber
BRAKES			
Type	Rear Drum	Rear Drum	Rear Drum
Parking Brake	Yes	Yes	Yes
BASKET			
Included	Yes	Yes	No
SPECIAL FEATURES			
Automotive style Tierod Steering	Yes	Yes	No

# Limited Lifetime Warranty

Knee Trekker Knee Walker frame is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser.

The lifetime warranty on the frame does not cover device failure due to owner misuse, negligence, or normal wear and tear and does not extend to non-durable components.

Brake components and hardware are warranted against manufacturer defect for one year. Wheels, rubber components, grips, cable, and the knee pad are warranted for 6 months as these parts are subject to normal wear and tear and will need periodic replacement.

If within such warranty period any such product shall be proven to be defective, the part will be repaired or replaced at Lighthouse's discretion.

Warranty covers full cost of replacement parts, it does not include labour or shipping costs.



If you have any questions  
call us at 1.866.609.6827





LightHouse  
Medical Ltd.



**1.866.609.6827**

[info@lighthousemedicaltd.com](mailto:info@lighthousemedicaltd.com)

[www.lighthousemedicaltd.com](http://www.lighthousemedicaltd.com)