

Anchor Bathing Chair

The Anchor provides a carefree bathing experience for children and caregivers alike. Safety and comfort are a priority so your child looks forward to bath time. We designed the anchor knowing your child's needs are unique and constantly evolving. That's why you can make angle and height adjustments with one hand on either side of the chair, all while your child is seated. It is also depth adjustable allowing for growth and long-term use. These agile and individualized functionalities ensure longevity, comfort, and hassle free bath time.



- Depth adjustable for growth, angle adjustable back, seat and calf support
- Height Adjustable
- One handed angle and height adjustments on either side of the chair
- Positioning straps and head support can be placed anywhere along the length of the chair
- Folds flat for easy storage and transport



ANCHOR Bath Chair		MSRP	Order
AN1200	ANCHOR Bath Chair XSmall - head support, positioning and leg belts included	\$900	
AN1400	ANCHOR Bath Chair Small - head support, positioning and leg belts included	\$950	
AN1600	ANCHOR Bath Chair Large - head support, positioning and leg belts included	\$1,000	
AN2100	ANCHOR Rolling Shower Base	\$950	
AN2200	ANCHOR Tub Stand	\$650	
AN2110	ANCHOR Lateral Support	\$99	
AN2111	ANCHOR Positioning Straps X-Large Pair	\$78	

	X-SMALL	SMALL	LARGE
SPECIFICATIONS			
Seat Width	12"	14.5"	16"
Seat Depth	10.5"-13"	10.5" - 13"	13"-16"
Back Height	18.5"	24"	32"
Back Width	14.5"	17"	18.5"
Calf Support Length	10.2"	11"	14"
Seat Height	5.75"-10.5"-14"-15"	5.75"- 11"-15.5"-17"	5.75"-13"-18.25"-20.5"
Back Angle	0°, 30°, 60°, 90°	0°, 30°, 60°, 90°	0°, 30°, 60°, 90°
Calf Support Angle	30°, 60°, 90°	0°, 30°, 60°	0°, 30°, 60°
Product Weight	12.75 lbs	15.5 lbs	18 lbs
Weight Capacity	155 lbs	155 lbs	155 lbs